

# ToRC Avidity Buffer

**Buffer Solution for Avidity testing of IgG class of antibodies for ToRC Diagnosis**  
FOR IN VITRO DIAGNOSTIC USE ONLY  
Store at 2°C to 8°C

## INTENDED USE

ToRC Avidity Buffer is intended for determining the percentage avidity binding of IgG class of antibodies in Toxoplasma, Rubella and CMV diagnosis using Bios Microwell ELISA reagent kits. For in Vitro Diagnostic Use only.

## INTRODUCTION

Avidity is a measure of antigen to antibody binding. Avidity Test helps in discriminating primary infection from secondary infection/re-infection when a sample reports positive for both IgG and IgM class of antibodies. Sometimes it is not sufficient to test only for IgM antibodies, as the presence of IgM class of antibodies may be due to the persistence of IgM antibodies from past infection or asymptomatic re-infection without risk for the foetus. For this reason, it is useful to assay the avidity of IgG antibodies. The low avidity of IgG antibodies is therefore an indication of recent or current infection.

## PRINCIPLE

ToRC Avidity Buffer acts as a denaturing agent and tends to dissociate the weaker IgG binding with antigen coated onto the microplate. However, if IgG antibodies are persistent against infection for long time a stronger bond is formed with antigen coated onto the microplate and doesn't get dissociated by Avidity buffer. This in turn helps in finding the strength of IgG antibodies in the patient's sample.

## STORAGE AND STABILITY


ToRC Avidity Buffer is stable at 2°C to 8°C up to expiry date printed on the label.

**NOTE:** ToRC Avidity Buffer is a ready to use reagent and to be used only with Bios Microwell ELISA reagents.

## TEST PROCEDURE FOR AVIDITY TESTING

- Place the desired number of coated strips into the holder. Consider 2 micro wells per sample along with NC and PC.
- Prepare 1:40 dilutions by adding 5 µl of the test samples, negative control & positive control to 200 µl of sample diluent and mix well.
- For the reagent blank, dispense 100 µl of sample diluent in A1 well position. Then dispense 100 µl of diluted negative control, positive control & serum samples in duplicates in appropriate wells, considering 1st well for Avidity Buffer and second well for wash Buffer. Tap the holder to remove air bubbles if any from the liquid. Incubate for 60 minutes at 37°C.
- Wash each well three times by filling approximately 350 µl diluted wash buffer & blot dry.
- Dispense 100 µl of Avidity Buffer in first well of each serum sample and controls and add 100 µl diluted Wash Buffer (1X) in the second well respectively and incubate for 10 minutes at 37°C.
- Wash each well three times by filling approximately 350 µl diluted wash buffer & blot dry.
- Dispense 100 µl of enzyme conjugate to each well and incubate for 45 minutes at 37°C.
- Wash each well three times by filling approximately 350 µl diluted wash buffer & blot dry.
- Dispense 100 µl of TMB Substrate to each well and incubate for 20 minutes at room temperature, away from direct light.

- Add 100 µl of Stop Solution to stop the reaction.
- Read O.D. at 450 - 630 nm with an ELISA reader.



## CALCULATION & INTERPRETATION OF THE RESULTS:

The percentage of avidity of the samples is expressed and calculated using a ratio between the O.D of wells containing the Avidity Buffer and those with Wash Buffer, subtracting the value of the blank.

$$\frac{\text{O.D with Avidity Buffer} - \text{O.D of Blank}}{\text{O.D with Wash Buffer} - \text{O.D of Blank}} \times 100 = \% \text{ of Avidity}$$

$$\text{O.D with Wash Buffer} - \text{O.D of Blank}$$

Ratio > 35% - High Avidity IgG  
Ratio < 30% - Low Avidity IgG  
Ratio 30-35% - Medium Avidity IgG

## IMPORTANT NOTE

- Inclusion of Blank with sample diluent is mandatory.
- This assay is a temperature sensitive assay. The best temperature condition for this assay is 37°C.

## BIBLIOGRAPHY

- Turune, H.J., P.O. Leinikke, and K. M. Saari. Demonstration of Intraocular Synthesis of Immunoglobulin G Toxoplasma Antibodies for

- Specific diagnosis of Toxoplasmic Chorioretinitis by Enzyme Immunoassay. J. Clin. Microbiol. 17:988-992, 1983.
- Lin, T.M., S.P. Halbert and G.R. O'Connor. Standardized Quantitative Enzyme-linked Immunoassay for Antibodies to Toxoplasma Gondii. J. Clin. Microbiol. Vol.11, 6:675-681, 1980.
- Gravell, M., P.H. Dorcett, O. Gutenson, and A.C. Ley. Detection of antibody to rubella virus by enzyme-linked immunosorbent assay. J. Infect. Dis. 136:S300, 1977.
- Hermann, K.L. Rubella virus. Manual of Clinical Microbiology, 3rd Edition. Lennette, Balows, Hausler, Truant (ed). Chapt. 86:862, 1980.
- Katz, S.L. Rubella (German measles). Zinssmer Microbiology, 18th Edition. Jolik, Willett, Amos (ed). Chapt. 75:1067, 1985.
- Voler, A., J.E. Bidwell, et al. Manual of clinical immunology. Chapter 69. Rose, N. and Friedman, H. eds. Am. Soc. Microbiol. p.506, 1985.
- Cremer, N.E. Antibodies in serodiagnosis of viral infection. p. 73. In Lennett E.H. ed. Laboratory diagnosis of viral infection. MerceL Dekker, Inc., New York, 1985.
- Starr, S.E. and H.M. Friedman. "Human CMV." Chapter 65. In Manual of Clin. Microbiol., 4th ed., Lennett, E.H. et al ed. Am. Soc. Microbiol. pp.771-719, 1985.
- Harry E Prince, Mary Lapé-Nixon; Focus Diagnostics Reference Laboratory, Cypress, California, USA. Role of cytomegalovirus (CMV) IgG avidity testing in diagnosing primary CMV infection during pregnancy.
- Department of Obstetrics and Gynecology, Kobe University Graduate School of Medicine, Kobe, Japan; Low IgG avidity and ultrasound fetal abnormality predict congenital cytomegalovirus infection.
- Murat Hökelek, MD, PhD Professor, Department of Clinical Microbiology, Istanbul University Cerrahpasa Medical Faculty, Turkey. What is the role of immunoglobulin G (IgG) avidity testing in the diagnosis of toxoplasmosis?

## SYMBOL KEYS

	Temperature Limitation		Consult Instructions for use
	Manufacturer		In vitro Diagnostic Medical Device
	Use by		Catalogue Number
	Date of Manufacture		Batch Number / Lot Number
	This side up		Contains sufficient for <-> tests
	Do not reuse		

Orchid Biomedical Systems (P) Ltd.\*

M 46-47, Phase III B, Verna Industrial Estate, Verna, Goa - 403 722, INDIA.  
Email id: sales@orchidbiomedical.com

\*A Revvity Inc. Company